

Systemstechnologies

The Experts in Wireless Nurse Call



Fall Prevention

Wired & Wireless Options for Bed Exit
& Fall Prevention

Quieter CordLess® – Wireless Fall Pad Monitor

Eliminates In-Room Noise

The FALTL-2100W wireless fall monitor can be mounted away from residents, eliminating in-room alarm noise. The monitor can be easily wall-mounted at the bedside or in a hallway up to 40 feet away from wireless fall pads using the provided hardware.

Patented Safety Features

This monitor features patented, cutting-edge technology that automatically lets a caregiver know if the signal between pad and monitor is lost or if the monitor is out of range. Featuring an optional patented TamperProof™ setting that, when enabled, requires caregivers use the provided Caregiver Key™ to reset the alarm, ensuring that residents cannot accidentally reset the monitor.

Wireless

This system has NO cords between pads and monitors, which reduces tripping hazards and eliminates broken or tangled pad cords. One monitor provides alerts for up to two wireless fall pads or one wired pad in the room. Available pad options are bed, chair, and floor.

Key Features

- Quiet. No sound to restrain residents.
- Wireless pads reduce tripping hazards and eliminates broken or tangled pad cords.
- Patented Pad Signal Lost notification alerts caregivers if sensor pad and monitor are no longer communicating.
- Optional TamperProof™ Caregiver Key allows control of who can reset the monitor.
- Raised visual status light enables caregivers to see alarm from a distance.
- Monitor up to two wireless pads at one time: Bed, chair, and floor.
- Three AA batteries included.
- Optional AC adapter available (AC-05).

Two Options for Bed Exit Equipment To Fulfill All Your Needs



Safety Auto-Reset Monitor (Wired Version)

SAFE w/ "Always ON" Technology

The FALTL-2100ST Monitor is designed for use with all Smart Caregiver corded sensor pad and floor mats. Simply plug a sensor pad into the monitor and place the pad under the resident in a chair or bed. When the resident gets up and pressure is removed, the alarm will sound, notifying caregivers that the resident is on the move. This alarm features a Safety Auto-Reset™ Function which will automatically reset and silence the alarm when pressure is re-applied to the sensor pad.

Optional TamperProof™ Setting

When set to TamperProof™, the monitor alert can ONLY be reset using the provided Caregiver Key™. Selecting this setting ensures that only the designated caregiver can silence the alarm.

Key Features

- Quiet - No sound to restrain residents.
- TamperProof™ Caregiver Key allows control of who is able to reset the monitor.
- Three AA batteries included.
- Optional AC adapter available (AC-02).
- Monitor different devices types: Wired bed, chair, or floor pad.
- Visual indicator.



Both wired and wireless options are compatible with our MicroVision and VisionLink II Nurse Call Systems.

Contact us today to find the best option for you!



Celebrating 30 years of life safety excellence!

Systemstechnologies

Toll-Free: 888-826-3394

Phone: 208-762-6800

Fax: 208-762-4877

11310 N. Government Way

Hayden, ID 83835

Website: www.wirelessnursecall.com

Email: sales@wirelessnursecall.com

